

IHE Integration Statement

Sectra RIS

Sectra RIS (Europe/Pacific), Version 22.2, September 2020



SECTRA

Knowledge and passion

Contents

1	Introduction	1
2	About IHE	2
2.1	Integrating the Healthcare Enterprise (IHE)	2
2.2	References	2
3	IHE integration statement	3

1 Introduction

The Integrating the Healthcare Enterprise (IHE) Integration Statement was developed as a high-level statement of integration for products, which were intended to interoperate, based on IHE Integration Profiles. This document identifies the IHE Actors, IHE Integration Profiles and options Sectra AB has implemented in each product version.

The Sectra IHE Integration Statement has been created to assist the reader in determining whether and to what extent interoperability with other vendors might be supported.

The IHE Technical Framework specifies a subset of the working mechanism of the healthcare enterprise and defines their exchanges in terms of a set of coordinated transactions. The actors and transactions expressed in the IHE Technical Framework are generalizations of actual healthcare information system environments. While some of the transactions are customarily carried out by specific product groups (e.g. HIS, RIS, PACS, or modalities), the IHE Technical Framework purposely avoids relating functions or actors with such product groups. For each actor, the IHE Technical Framework identifies only those functions linked with integrating information systems. The IHE definition of an actor should not be interpreted as the absolute definition of any product that might implement that particular actor. In addition, the framework itself should not be taken as the absolute definition of healthcare information system architecture.

This IHE Integration Statement affords the reader with a sophisticated view of supported IHE integration profiles. For further exploration, additional information can be found in the DICOM and HL7 Conformance Statements of this product and in the IHE Technical Framework.

This IHE Integration Statement does not guarantee successful interoperability of this product with other vendor's products. It is the user's responsibility to thoroughly analyze the application requirements and to specify a solution that integrates the Sectra PACS or Sectra VNA with other vendor's software and systems.

2 About IHE

The following topics are included in this chapter:

- [Integrating the Healthcare Enterprise \(IHE\)](#)
- [References](#)

2.1 Integrating the Healthcare Enterprise (IHE)

IHE is a multi-year initiative undertaken by medical specialists, administrators, information technology professionals and manufacturers, sponsored by the Radiological Society of North America (RSNA) and the Healthcare Information and Management Systems Society (HIMSS), whose sole purpose is to improve the way computer systems in healthcare share information.

IHE is not a standards organization. Instead it promotes coordinated use of existing standards such as DICOM and HL7 to develop workflow solutions for the healthcare enterprise. Systems designed in agreement with IHE profiles communicate with one another better, are easier to implement, and facilitate the efficient access of information. Physicians, nurses, administrators and other healthcare professionals foresee a day when vital information can flow seamlessly from system to system and be readily available at the point of care. IHE is intended to make this vision a reality by improving the condition of systems integration and eliminating barriers to optimal patient care.

Though begun in the USA, the IHE is a world wide initiative. The Regional IHE Initiatives currently encompass IHE – Asia Oceania, IHE – EU (Europe) and IHE – NA (North America). The National and Regional IHE Initiatives coordinate the deployment of the IHE Technical Framework and Integration Profiles and give input to their development related to local needs and issues. They are sponsored and overseen by relevant national/regional professional associations and include membership and staff of these groups as well as local vendor representatives. These groups may rely on funding from other organizations (e.g., governmental bodies or foundations) and are responsible for securing and managing any such funds.

2.2 References

General IHE information:

- IHE home: <http://www.ihe.net>

Vendor's IHE information:

- <https://sectra.com/medical/knowledge-center/conformance-statements/>

3 IHE integration statement

Vendor	Product Name	Version	Date
Sectra AB	Sectra RIS (Europe/Pacific)	22.2	2020-09-22

This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:

Integration Profile Implemented	Actors Implemented	Options Implemented
Scheduled Workflow	Department System Scheduler/Order Filler	None
	Performed Procedure Step Manager	None
Patient Information Reconciliation	Department System Scheduler/Order Filler	None
	Performed Procedure Step Manager	None
Presentation of Grouped Procedures	Department System Scheduler/Order Filler	None
	Performed Procedure Step Manager	None
Consistent Time	Time Client	None
Audit Trail and Node Authentication	Secure Application	None
Cardiac Catheterization Workflow	Department System Scheduler/Order Filler	None
	Performed Procedure Step Manager	None
Echocardiography Workflow	Department System Scheduler/Order Filler	None
	Performed Procedure Step Manager	None


LABEL

Product Sectra RIS (Europe/Pacific)

Version 22.2

Manufacturer Sectra AB
Teknikringen 20
SE-58330 Linköping
Sweden
www.sectra.com

Contact <https://sectra.com/medical/about-sectra/contact/>



Regulatory Clearance Statement

The quality system of Sectra AB [Sectra] conforms to ISO 9001, ISO 13485 and ISO 27001. All Sectra medical devices have obtained regulatory clearance for those markets where Sectra sells and deploys its devices, e.g. EEA, USA, Canada, Australia. For further regulatory information, please contact Sectra.

Disclaimer

Sectra AB is not responsible for problems caused by changes in the operating characteristics of the computer hardware or operating system, which are made after the delivery of the software, or for problems that occur as a result of the use of Sectra software in conjunction with non-Sectra software other than the software explicitly covered in this documentation. For a complete description of warranty, refer to the End User License Agreement supplied with the Sectra product.

Trademark and Patent Disclaimer

Sectra, Sectra Imtec and the Sectra logotype, are registered trademarks of Sectra AB. FIMAG and Image Exchange Portal are registered trademarks of Burnbank Systems Ltd in the UK. IDS7, Sectra BizTrack, Sectra CloudFlex, Sectra DoseTrack, Sectra IEP, Sectra Image Central, Sectra Image Lab, Sectra LiteView, Sectra UniView, Sectra One Connect, Sectra Open Archive, Sectra VNA, Sectra UserInfluence, and WISE are trademarks of Sectra Imaging IT Solutions AB. Sectra Preop Online and Sectra OneScreen are trademarks of Sectra AB. RxEye is a registered trademark of RxEye AB.

Windows is a registered trademark of Microsoft Corporation in the United States and other countries. *Mac*, *Safari*, and *iPad* are registered trademarks of Apple Inc. in the United States and other countries. All other names/products by ® or ™ are registered trademarks of the respective manufacturer.

The intellectual property of Sectra includes a number of patents, for full information refer to: <https://sectra.com/medical/about-sectra/patents/>

Copyright and Legal Statement

© Sectra AB, Sweden, 2020

All rights are reserved. Reproduction or transmission in whole or part, in any form or by any means, electronic, mechanical or otherwise, is prohibited without written consent of the copyright owner.

Copyrights and all other proprietary rights in any software and related documentation ("Software") made available to you rest exclusively with Sectra AB. No title or ownership in the Software is conferred to you. Use of the Software is subject to the end user license conditions as are available on request.

To the maximum extent permitted by law, you shall not decompile and/or reverse engineer the software or any part thereof.

Sectra AB
Teknikringen 20
SE-58330 Linköping
Sweden
Phone: +46 13 23 52 00
E-mail: info.medical@sectra.com

SECTRA

Knowledge and passion