


Echo Viewer

Sectra's Echo Viewer transforms the viewing of ultrasound series, delivering unparalleled speed and seamless interaction. It is available both in Sectra IDS7 and in the zero-footprint web viewer Sectra UniView. Built for easy-to-use echocardiography reviewing; it loads exams instantly, making it easy to navigate and ensuring a lag-free playback to save time.



Features

- Made for Ultrasound viewing
- Smooth scrolling experience
- Displaying 16 series simultaneously
- High frames-per-second rates
- Comparison mode
- Viewing of Angiography images
- Viewing of IVUS images

Functionalities

- Distance measurement
- Area measurement
- Basic doppler measurements

Coming soon!

- TTE standard measurements
- Cardiovascular measurements
- Echo workflow supported by AI
- Import, export, and display DICOM SR measurements

Read about Sectra Cardiology Imaging:
medical.sectra.com/cardiology